

PHOTOSYNTHETIC PIGMENT DETERMINATION USING ETHANOL EXTRACTION

Principle

A water sample of known volume containing photosynthetic organisms is filtered through glass fibre filter paper and chlorophyll is extracted in a refrigerator overnight or at 60°C for 1 hour in ethanol (preferred over methanol for safety reasons). The absorbance of the chlorophyll extract is read spectrophotometrically and readings used to calculate the concentration of chlorophyll *a*, phaeophytin and carotenoid accessory pigments.

Procedures

1. Filter 1 litre of water (volume depends on the amount of suspended material) through a pre-weighed Whatman 47 mm GF/C filter paper, using slight vacuum (no more than 0.5 Bar=30cm. Hg).
2. Weigh filter and place in a stoppered centrifuge tube, eg Hach tube. Freeze filter, or for preference put immediately into extractant. (Water content of the filter determined by weighing the filter dry and after filtration: normally is ~0.4-0.5ml)
3. Extract each filter in 19 volumes of 100% ethanol (industrial methylated spirit); eg. if water content was 0.5g, add 9.5 ml. (This procedure ensures that the water content of the final mix is not greater than 5%). Mark the solvent level on the tube, invert several times to mix.

Solutions : 100% and 95% ethanol (industrial methylated spirit) saturated with MgCO₃

4. Wrap tubes in aluminium foil to exclude light: light degrades chlorophyll.
5. Place tubes in a waterbath at 60°C for 1 hour. This speeds extraction and destroys chlorophyllase.
6. Remove tubes, invert several times, and allow to cool on ice. They can be stored overnight in a fridge or freezer. Any evaporative loss should be made up with 100% ethanol, topping up to the mark again. (extraction for 24 hours is recommended)
7. Centrifuge without removing the filter at full speed (about 3500 rpm) in a bench centrifuge for 15 minutes. Decant or pipette out the extract: check volume. Alternatively the filter may be removed (squeeze out the solvent on the side of the tube) before spinning. The volume of solvent must be greater than 8 ml to read in a 4 cm path length spectrophotometer cell. For convenience the volume may be adjusted to exactly 10 ml; otherwise record the volume to the nearest 0.1 ml for use in subsequent calculation.

Cells used should be read against the reference, containing the same blank solution*, before use at all wavelengths, as not all cells are perfectly matched: subtract any differences from the extract readings. If acidifying extracts, read acidified blank cells also.

8. Read the absorbance against a suitable reference [*MgCO₃-saturated 95% ethanol which has been centrifuged or filtered (Whatman No. 1)] at the following wavelengths, in a 4 cm glass cell : 750, 665, 510 and 480 nm. Return wavelength to 665.
The OD at 665 should be between 0.2 and 0.8: it may be necessary to dilute more concentrated samples (or read them in a shorter path length cell) -see Marker et al, 1980.
9. To determine phaeophytin **
Acidification of chlorophyll degrades it to phaeophytin with a change in Specific Absorbance Coefficient (SAC), which can be used to estimate the amount of phaeophytin present before acidification. In methanol (ethanol) the extract must be neutralised before re-reading absorbance, in order to prevent complications due to differences in pH. (pH after neutralisation can be checked with pH paper)

The method assumes an extract volume of 10 ml.

- (a) Add directly to the cuvette 0.1 ml of 1 M HCl (10%), mix.
- (b) Leave for 2 minutes
- (c) Add 0.13 ml* of 1 M 2-phenylethylamine in ethanol (12.1g in 100 ml), mix.
- (d) Read absorbance at 665 and 750 nm.

* with new solutions, check the amount of base needed to neutralise 0.1ml acid: add 0.1ml base and check by spotting on pH paper, then add 10µl aliquots until pH7 is reached

10. Calculations:

Subtract cuvette blanks at each wavelength from all readings.

Subtract A_{750} (turbidity blank) from all other readings.

Equations refer to corrected values: (Parsons, Maita and Lalli (1984) equations for acetone^{**})

$$\text{chlorophyll } a \text{ } (\mu\text{g l}^{-1}) = \frac{29.11^{**} [A_{665_o} - A_{665_a}] s d}{V \rho}$$

^{**} modification to factor if used in ethanol: change to 29.11 from 26.7.
Not rigorously tested yet

$$\text{phaeophytin } a \text{ } (\mu\text{g l}^{-1}) = \frac{29.11^{**} [1.7 (A_{665_a}) - A_{665_o}] s d}{V \rho}$$

d = dilution factor (if necessary to reduce A665 to <0.8) eg. for 1 in 5 dilution, d = 5

ρ = path length of cell, cm

s = solvent extract volume, ml

V = sample volume, litres

A_{665_o} = absorbance at 665 nm before acidification

A_{665_a} = absorbance at 665 nm after acidification and neutralisation.

Notes

Method for phaeophytin is according to the recommendations of Marker, Crowther and Gunn (1980) and Marker *et al.* (1980) for methanolic extracts. Webb *et al.* (1992) found that spectrophotometric methods greatly overestimate phaeophytin compared to chromatographic methods, and therefore do not recommend its use. As an alternative, calculate chlorophyll a according to this equation:

$$\text{Chlorophyll } a \text{ } (\mu\text{g l}^{-1}) = \frac{11.99 (A_{665} - A_{750}) s}{V \rho}$$

Marker *et al.* (1980)
Recommended equation for ethanol
(cf. constant in methanol=12.99)

Carotenoids: calculate as 'microscopic pigment units' = μspu .

Most formulae given refer to acetone extracts, e.g. Strickland and Parsons (1968):

(a) $\mu\text{spu l}^{-1} = \frac{4.0 (A_{480} - A_{750}) s}{V \rho}$ if predominantly Chlorophyta or Cyanophyta

(b) $= \frac{10.0 (A_{480} - A_{750}) s}{V \rho}$ if predominantly Chrysophyta or Pyrrophyta

Parsons, Maita and Lalli (1984) give this formula:

(c) $\mu\text{spu l}^{-1} = \frac{7.6 [(A_{480} - A_{750}) - 1.49 (A_{510} - A_{750})] s}{V \rho}$

[Foy (1987) used formula (a) with methanolic extracts of Cyanophyta.]

NB. Extraction from stones:

Drain excess moisture from the stones and place each stone in a resealable plastic bag (double-wrap for security), weigh.

Add an appropriate volume of extractant for the size of the stone, seal and re-weigh. Avoid exposure to bright light. Appropriate volumes:- eg. 50 ml., depending on size of stone and amount of algal cover. With dense mats estimate the water content and use 9.5 volumes of 100% alcohol, diluting further with 95% if needed.

Place in a water bath at 60°C for 1 hour, cool (on ice). Re-weigh to check no evaporative loss (or gain though leakage).

Remove 10ml to a Hach tube, centrifuge. Extracts may require dilution: A_{663} should not exceed 0.8.

Measure absorbance as with filtered samples.

Calculate chlorophyll yield per unit area (cm²) of stone surface. Stone surface covered by algae can be determined by wrapping the area with foil and comparing the foil weight with a known area. Area takes the place of sample volume (V) in the formula.

REFERENCES

- Arvola, L. 1981. Spectrophotometric determination of chlorophyll-a and phaeopigments in ethanol extractions. *Ann Bot Fennici*. **18**. 221-227
- Foy, R H. 1987. A comparison of chlorophyll a and carotenoid concentrations as indicators of algal volume. *Freshwater Biology*. **17(2)**. 237-250
- Hansson, L-A. 1988. Chlorophyll-a determination of periphyton on sediments: identification of problems and recommendation of method. *Freshwater Biology*. **20**. 347-352
- Holm-Hansen, O and Riemann, B. 1978. Chlorophyll a determination: improvements in methodology. *Oikos*. **30**. 438-447
- Jespersen, A-M and Christoffersen, K. 1987. Measurements of chlorophyll-a from phytoplankton using ethanol as extraction solvent. *Archives of Hydrobiology*. **109(3)**. 445-454
- Lorenzen, C J. 1967. Determination of chlorophyll and pheo-pigments: spectrophotometric equations. *Limnology and Oceanography*. **12**. 343-346
- Marker, A F H. 1972. The use of acetone and methanol in the estimation of chlorophyll in the presence of phaeophytin. *Freshwater Biology*. **2**. 361-385
- Marker, A F H. 1977. Some problems arising from the estimation of chlorophyll a and pheophytin a in methanol. *Limnology and Oceanography*. **22(3)**. 578-579
- Marker, A F H, Crowther, C A and Gunn, R J M. 1980. Methanol and acetone as solvents for estimating chlorophyll-a and phaeopigments by spectrophotometry. *Archives of Hydrobiology Bulletin (Ergebnisse der Limnologie)*. **14**. 52-69
- Marker, A F H, Nusch, E A, Rai, H and Riemann, B. 1980. The measurement of photosynthetic pigments in freshwaters and standardization of methods: conclusions and recommendations. *Archives of Hydrobiology Bulletin (Ergebnisse der Limnologie)*. **14**. 91-106
- Moss, B. 1967 (a). A spectrophotometric method for the estimation of percentage degradation of chlorophylls to pheo-pigments in extracts of algae. *Limnology and Oceanography*. **12**. 335-340
- Moss, B. 1967 (b). A note on the estimation of chlorophyll-a in freshwater algal communities. *Limnology and Oceanography*. **12**. 340-342
- Nusch, EA. 1980. Comparison of different methods for chlorophyll and pheopigment determination. *Archives of Hydrobiology Bulletin (Ergebnisse der Limnologie)*. **14**. 14-36
- Parsons, T R and Strickland, J D. 1963. Discussion of spectrophotometric determination of marine plant pigments with revised equations for ascertaining chlorophylls. *Journal of Marine Research*. **21**. 155-163
- Parsons, TR, Maita, Y & Lalli, CM. 1984. *A Manual of Chemical and Biological methods for seawater analysis*. Pergamon Press, Oxford.
- Rai, H. 1980. Comparison between trichromatic spectrophotometric equations: sources of error. *Archives of Hydrobiology*. **88(4)**. 514-517
- Strickland, J D H *et al.* 1968. Spectrophotometric determination of chlorophylls and total carotenoids. *Fisheries Research Board of Canada Bulletin*. **167**. 185-192
- Strickland, J D H and Parsons, T R. 1972. *A Practical Handbook of Seawater Analysis*, 2nd Edition. Fisheries Research Board of Canada, Ottawa. .
- Tett, P, Kelly, M G and Hornberger, G M. 1975. A method for the spectrophotometric measurement of chlorophyll a and pheophytin a in benthic microalgae. *Limnology and Oceanography*. **20**. 887-896
- Tett, P, Kelly, M G and Hornberger, G M. 1977. Estimation of chlorophyll a and pheophytin a in methanol. *Limnology and Oceanography*. **22**.579-580
- Webb, D J, Burnison, B K, Trimbee, A M and Prepas, E E. 1992. Comparison of chlorophyll-a extractions with ethanol and dimethyl sulfoxide-acetone, and a concern about spectrophotometric phaeopigment correction. *Canadian Journal of Fisheries and Aquatic Sciences*. **49(11)**. 2331-2336
- Wood, L W. 1985. Chloroform-methanol extraction of chlorophyll a. *Canadian Journal of Fisheries and Aquatic Sciences*. **42**.38-43