

## Total Nitrogen analysis

Total nitrogen (or total soluble, if the sample is filtered) can be assayed by oxidation with alkaline persulphate to nitrate, followed by analysis using the hydrazine reduction technique [modified from APHA, 1995].

Reagents:

### Digestion:

#### Digestion reagent:

Potassium persulphate  $K_2S_2O_8$  (AR) 20.1 g  
Sodium Hydroxide NaOH (AR) 3.0 g  
Deionized water to 1000 ml  
Make appropriate volume just prior to use

#### Borate buffer:

Boric acid  $H_3BO_3$  (AR) 61.8 g  
Sodium Hydroxide NaOH (AR) 8.0 g  
Deionized water to 1000 ml

### Digestion check Standards:

#### Glutamic acid stock solution:

dry glutamic acid at 105°C for 24h

Dissolve 1.051 g in deionized  $H_2O$ , dilute to 1000 ml 1.00 ml = 100 µg N  
preserve with 2 ml  $CHCl_3$

Glutamic acid intermediate solution: Dilute stock 1 in 10 to give 10 µg per ml  
for digestion check standard, dilute 29 ml of stock to 100 ml = 2.9 µg.ml<sup>-1</sup>

#### [Yeast or Yeast Extract:

Suspend/dissolve 102 mg in 100 ml  $H_2O$ : total N should be 2 mg.l<sup>-1</sup>. [ total P in this standard should be 3.97 mg.l<sup>-1</sup>. ] This has been tested as a digestion check standard, but the results are equivocal]

### Nitrate assay:

NaOH reagent: 250 ml *caution: caustic*

NaOH 25 g, add to 237.5 ml deionized  $H_2O$ , store in air tight plastic bottle at room temp. for 1 month.

Reducing reagent: *caution: poisonous, highly reactive*

A: Hydrazine sulphate 2.7 g in 100 ml deionized  $H_2O$  *caution: poisonous*  
(heat to dissolve)  
B: Copper sulphate ( $CuSO_4 \cdot 5H_2O$ ) 0.25 g in 100 ml deionized  $H_2O$  *caution: harmful*  
C: Zinc sulphate ( $ZnSO_4 \cdot 7H_2O$ ) 5.3 g in 100 ml deionized  $H_2O$  *caution: harmful*

Store all solutions at room temp: stable indefinitely until mixed.

On day of use mix: 27.4 ml  $H_2O$ ; 10 ml A; 1.5 ml B; 1.1 ml C.  
(do not use if more than 1 day old)

Colour Reagent: add in the following order:

deionized  $H_2O$  90 ml  
Orthophosphoric acid 50 ml  
Sulphanilamide 5 g  
N-1-naphthylethylenediamine dihydrochloride (NEDD) 0.25 g  
dissolve and add (use for rinsing reagents in, etc.) deionized  $H_2O$  57.5 ml  
store refrigerated in a dark bottle: renew monthly

### Standard Nitrate solution:

KNO<sub>3</sub> (anhydrous) 7.221g

Dissolve in deionized H<sub>2</sub>O to 1000ml (volumetric flask) @ 20°C: add 2 ml Chloroform as preservative: keep for no longer than 1 week.

This stock is 1000mg/l NO<sub>3</sub> - N. Dilute to 5000µg.l<sup>-1</sup> for working strength: this solution should not be stored. make up standards for analysis up to 3000µg.l<sup>-1</sup>.

### **Method:**

- 1 add 20 ml of standard solutions or sample (suitably diluted) to pre-weighed 50-100 ml conical flasks (cleaned first by autoclaving with digestion reagent: thereafter by acid-washing)
- 2 add 10 ml of digestion solution, mix. Cover flasks with small glass beakers.
- 3 Place in autoclave: heat for 30 minutes at 121 °C.
- 4 Allow to cool to room temperature, check volume (by weighing) and make up to 30 ml if necessary: add 2 ml borate buffer, mix.
- 5 Distribute 5 ml aliquots from each flask to 3 acid-washed glass test tubes .
- 6 Add 725µl NaOH reagent, mix (vortex).
- 7 Add 420µl of freshly made reducing agent and mix.
- 8 Wait 4-5 minutes (or more)
- 9 Add 725µl Colour reagent, mix.
- 10 Leave for 30 min for full colour development, mix and read at 535nm.  
Absorbance of a 1000µg/l standard in a 1 cm cell should be approx. 0.5-0.7.

### References

APHA. 1995. Standard Methods for the Examination of Water and Wastewater, 19<sup>th</sup> edition. American Public Health Association.

D'Elia, CF and Steudler, PA. 1977. Determination of total nitrogen in aqueous samples using persulfate digestion. Limnology and Oceanography. 22. 760-764.

Dilutions for standard curve: volumes to make 20ml for digestion reaction

Concentration µg/l	5 mg/l stock	Water
0	-	20ml
500	2ml	18
1000	4	16
1500	6	14
2000	8	12
2500	10	10
3000	12	8

## **Total Phosphorus Analysis:**

by acid persulphate digestion followed by ascorbic acid procedure.

For Total Soluble P (TSP), filter the water sample before analysis, otherwise the procedure gives Total P (TP)

### **Reagents:**

30% H<sub>2</sub>SO<sub>4</sub>                      30 ml conc H<sub>2</sub>SO<sub>4</sub> added to 50ml deionised water, made up to 100 ml when cool.

*Caution highly corrosive, great heat emitted during dilution: dilute into beaker in a basin of cold water, in fume hood.*

8% potassium persulphate:              8g potassium persulphate dissolved in 100 ml deionised H<sub>2</sub>O, stirring on a hotplate to dissolve. Prepare just before use and dispense while still warm.

### **Procedure:**

- 1            Add 5 ml of sample, standards, and blank to re washed Hach COD tubes, in triplicate.
- 2            Add 0.1 ml 30% H<sub>2</sub>SO<sub>4</sub> and 0.5 ml fresh potassium persulphate solution.
- 3            Put screw caps on tubes and tighten
- 4            Place in autoclave and autoclave for 30 min at 15 psi (121°C)
- 5            Allow to cool
- 6            Add 1 ml of freshly prepared mixed reagent (see Ascorbic acid method) and vortex to mix.
- 7            Leave 20-30 min, read absorbance at 882nm in a 1cm flowcell zeroed against a deionised water blank